

SUMMIT LASER SPECIFICATIONS SHEET



7 Watt

10 Watt

15 Watt

20 Watt

| | |
|------------------------|--|
| Description | Desktop – Summit Series Class IV Laser |
| Model | Summit Series 7 Watt / 10 Watt / 15 Watt / 18 Watt |
| SKU | AL-G7x / AL-G10x / AL-G15x / AL-G18x |
| Laser Type | GaAlAs Diode Laser |
| Wavelength, Laser | 980nm |
| Maximum Power | 7 watts / 10 watts / 15 watts / 18 watts |
| Operation Modes | Continuous Wave (CW) and Pulse |
| Pulse Duration | 10µs-3s |
| Repetition Rate | 0.2Hz-20KHz |
| Fiber Sizes | 550µM (Therapy) |
| Fiber Connector | SMA905 Connector |
| Pilot Beam | Red Diode Laser of 635nm, Power<5mW |
| Control Mode | True Color Touch Screen<5mW |
| Voltage/Current Rating | 110/220 VAC, 5A , 50/60 Hz |
| Dimensions | 215(W) x 245(L) x 315(H) mm |
| Weight | 9 lb. |
| Safety Compliance | CE0197 |

Specifications subject to change without notice.

